

Lounge Chair Series

Models:

BD7020, GA1954, GA1954X, SA1403, SA1403X, MB1954X, MB1403X

Assembly/Usage Instructions

Install Legs- (excluding MB1954X & MB1403X)

- Place the 1 ¼" washers on the ¼" x 3" bolts and insert through pre-drilled leg holes on sides and ends of the laminated base.
- Place steel leg on inside of base and over bolts. Place ¼" washers and lock nuts on to bolts and tighten with a 7/16" wrench. Be careful not to over tighten nuts.

Install Arms-

- Remove locking nut and small washer from back of chair, place arm bracket on to studs leaving large washer between arm and chair, replace washer and lock nut (hand tighten only at this time).
- Place arm post into arm assembly - level armrests and tighten arm bracket fully with a ½" wrench. Adjust armrest to desired level and secure by tightening arm adjuster handle.

Attach Lift Bar- (1954 & 1403 Lounge Series Only)

- Secure "U" shaped lift handle to foot of chair using 1/4" x 3 1/2" bolts with flat-head screwdriver.

Usage Instructions

Arm Adjustment-

The armrest is specifically designed to lift up and out of the way to allow easy entrance and exit by patient.

- To adjust the height of the armrests loosen the black adjuster handle and move the arm up/down or left/right. Re-tighten handle to secure arm into position. The black adjuster handle position can be adjusted by pressing in gray handle button and rotating handle.

Recline-

- **When bed is equipped with a hydraulic lift mechanism-** Place hands on lift handle and squeeze black release handle. The lounge can be stopped and held in any position by releasing the release handle at the desired position.
- **When bed is not equipped with a hydraulic lift mechanism-** Place hands on end of bed and slowly lift up.